

Blood type ID Card Centrifuge I IDC-3000



Blood type ID Card Centrifuge IIDC-3000 with standard programs are built-in for blood type determination and all kinds experiment of serology, which standardize the test and research of clinical. Users can set the programs with simple touch screen button. Suitable for blood group testing, micro column gel immune detection. Spins sample volumes up to 24 cards micro gel rotor with maximum speed 3000 rpm. Driven by DC brushless motor. Auto-balancing centrifuge for blood bank is particularly designed for the quick centrifuge in blood bank with storages up to 9 programs records for users. For safety we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function.

Accessories



Technical Specification:

Model	IIDC-3000
Speed Range (rpm)	3000 rpm
Max RCF	1450xg
Capacity	24 cards
Timer	1 to 99 min
Net Weight (kg)	16 kg
Electrical Requirements	AC220v/50Hz
Catalogue No.	1402150502

Technical Data:

Model	IIDC-3000
Capacity	24 cards
Max RCF	1450xg
Speed Range	3000 rpm
Noise	? 60dB?A?
Timer	1 to 99 min
Weight	16 kg
Electrical Requirement	AC220v/50Hz
Catalogue No.	1402150502



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States