

Chest Type Medical Freezers

-25°C ICMF-248



Chest Type Medical Freezers -25°C ICMF-248 design for a large amount of freezer storage, the available maximum capacities are 485 L with a temperature range of -10°C to -25°C by offering excellent temperature uniformity and stability. Build-in inner door to avoid more cooling leaking out, and provide more stable inside temperature performance and complete temperature alarm system featuring. Used to stores vaccines, blood plasma, and many other biological products. Built-in wheel casters to provide adjustable support safest and secure unit on any rough surface. Features cooling fan for instrument safety, durability, and endurance purpose. The range of medical freezers storage products will permit you to choose the perfect match to meet your requirements through cost and energy saving.

Technical Specification:

Model	ICMF-248
Capacity (liters)	485 L
Cabinet Type (Upright/Chest)	CHEST
Temperature Range (°C)	-10 to -25°C
Refrigeration System	HC-R290
Defrost	Manual
Alarms	High/ low temperature, Power off, Sensor malfunction, Temperature controller malfunction
Material	Inner : Stainless steel; Outer: Color sprayed steel
Overall Dimensions (WxDxH) cm	1655×755×915 mm
Certifications/Compliance	ISO9001□13485□14001,CE
Net Weight (kg)	124 kg
Gross Weight (kg)	138 kg
Electrical Requirements	220V/50Hz
Catalogue Number	500301002

Technical Data:

Model	ICMF-210	ICMF-248
Capacity	105 L	485 L
Temperature Range	-10 to -25?	
Materials	Inner : Stainless steel; Outer: Color sprayed steel	
Alarms	High/ low temperature, Power off, Sensor malfunction, Temperature controller malfunction	
Cooling	Direct cooling	
Class	CHEST	
Condenser	Wire-and-tube	
Sensor	NTC	
Refrigerant	HC-R290	
Evaporator	Submerged coil evaporator	
Defrosting	Manual	
Casters	4	
Compressor	CUBIGELN1	
Temperature Controller	Microprocessor	
Ambient Temperature	10 to 32?	
Certification	"ISO9001?13485?14001,CE"	
Power Consumption	120W	295W
Internal Dimension	480×380×665 mm	1455×525×665 mm
External Dimension	680×580×935 mm	1655×755×915 mm
Net Weight	62 kg	124 kg
Gross Weight	70 kg	138 kg
Electrical Requirement	220V/50Hz	
Catalogue No.	500300502	500301002



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States