

Upright Deep Freezers -40°C

IUDF-232



Upright Deep Freezers -40°C IUDF-232 is well organized single door low-temperature freezer as they contain shelves and draws, letting foods to be quickly sorted. The maximum capacities up to 322 L with a temperature range of -30°C to -40°C.

Equipped with automatic magnetic sealed door closes to ensure and heightens the refrigeration performance and operational reliability, and password protection to control operating parameters from being modified, robust and durable framing, antibacterial powdered inner and outer surface, refrigeration system uses highly efficient and eco-friendly and fluorine-free refrigerant. Furnished with a rechargeable battery, which can display the main information such as alarm and the internal temperature of the cabinet when the main power is cut off Built-in mechanical lock to ensure safety

Technical Specification:

Model	IUDF-232
Capacity (liters)	322 L
Cabinet Type (Upright/Chest)	Upright
Temperature Range (°C)	-30 to -40°C
Temperature Accuracy (°C)	0.1°C
Alarms	High/low temperature alarm, sensor error alarm, power failure, door opening
Material	Outer: White plastic coated galvanized steel; inner: stainless steel
Refrigerant Type	Non HCFC
Shelves	3
Overall Dimensions (WxDxH) cm	673 x 676 x 1886 mm
Casters	4
Net Weight (kg)	113 kg
Electrical Requirements	AC220V/50Hz
Catalogue Number	500351002

Technical Data:

Model	IUDF-260	IUDF-232
Capacity	263 L	322 L
Temperature Range	-30 to -40°C	
Materials	Outer: White plastic coated galvanized steel; inner: stainless steel	
Alarms	High/low temperature alarm, sensor error alarm, power failure, door opening	
Type	Upright	
Display	Digital	
Refrigerant	Non HCFC	
Shelves	3	4
Casters	4	
Interface		
Accuracy	0.1°C	
Certification		
Internal Dimension	508 x 455 x 1137 mm	508 x 455 x 1393 mm
External Dimension	673 x 676 x 1630 mm	673 x 676 x 1886 mm
Weight	98.5 kg	113 kg
Electrical Requirement	AC220V/50Hz	
Catalogue No.	500350502	500351002



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States