

Low Temperature Chest Freezers

-60°C ILCF-2105C



Low-Temperature Chest Freezers -60°C ILCF-2105C design for large amounts of freezer storage, the available maximum capacities: 105 L with a temperature range of -25°C to -60°C by offering excellent temperature uniformity and stability. Build-in inner door to avoid more cooling leaking out, and provide more stable inside temperature performance and complete temperature alarm system featuring. Used to stores vaccines, blood plasma, and many other biological products. Features cooling fan for instrument safety, durability, and endurance purpose. The range of medical freezers storage products will permit you to choose the perfect match to meet your requirements through cost and energy saving.

Technical Specification:

Model	ILCF-2105C
Capacity (liters)	105 L
Cabinet Type (Upright/Chest)	Single door chest
Temperature Range (°C)	-25 to 60°C
Refrigeration System	CFC-Free
Defrost	Manual
Alarms	High/low temperature, Sensor error, Thermostat failure
Material	Inner : Stainless steel; Outer: Color sprayed steel
Condenser	Copper tube and finned
Overall Dimensions (WxDxH) cm	680×580×935 mm
Certifications/Compliance	ISO9001□13485□14001,CE
Net Weight (kg)	62 kg
Gross Weight (kg)	70 kg
Electrical Requirements	220V/50Hz
Catalogue Number	500600502

Technical Data:

Model	ILCF-2105C	ILCF-2485D
Capacity	105 L	485 L
Temperature Range	-25 to 60?	
Materials	Inner : Stainless steel; Outer: Color sprayed steel	
Alarms	High/low temperature, Sensor error, Thermostat failure	
Cooling	Direct cooling	
Class	Single door chest	Double door chest
Condenser	Copper tube and finned	
Sensor	PTC	
Refrigerator	CFC-Free	
Evaporator	Submerged coil evaporator	
Defrosting	Manual	
Casters	4 to 6	
Compressor	CUBIGEL/1	
Temperature Controller	Microprocessor	
Ambient Temperature	10 to 32?	
Certification	"ISO9001?13485?14001,CE"	
Power Consumption	280 W	450 W
Internal Dimension	480×380×665 mm	1455×525×665 mm
External Dimension	680×580×935 mm	1655×755×915 mm
Net Weight	62 kg	124 kg
Gross Weight	70 kg	138.5 kg
Electrical Requirement	220V/50Hz	
Catalogue No.	500600502	500601002



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States