

Double Door Medical Refrigerator IUMR-2700D



Double Door Medical Refrigerator 2 to 20°C IUMR-2700D with dual-layer opening and 10 shelves for easy observation and storage of samples. The IUMR-2000D comes with the same design and features as IUMR-1000 having double doors instead.

The maximum capacity of 997 L and defrost automatically. Build-in inner fan-driven and forced-draft air circulation system provides excellent temperature uniformity. Build-in microprocessor controller to certify tolerances of internal temperatures. Equipped with buzzer and flashing lights, power failure alarm, doors open alarm, high and low-temperature alarm.

Technical Specification:

Model	IUMR-2700D
Capacity (liters)	997 L
Cabinet Type (Upright/Chest)	Double door Upright
Temperature Range (°C)	2 to 20°C
Refrigeration System	Non HCFC
Alarms	Power failure alarm, Temperature sensor failure alarm, High / Low temperature alarm, Door open alarm
Material	White plastic coated galvanised steel
Shelves	10
Overall Dimensions (WxDxH) cm	1220 x 872 x 1885 mm
Net Weight (kg)	207 kg
Electrical Requirements	AC220 V/ 50Hz
Catalogue Number	500301002

Technical Data:

Model	IUMR-2656D	IUMR-2700D
Capacity	662 L	997 L
Cabinet Type	Double door Upright	
Temperature Resolution	0.1°C	
Temperature Range	2 to 20°C	
Alarms	Power failure alarm, Temperature sensor failure alarm, High / Low temperature alarm, Door open alarm	
Refrigerant	Non HCFC	
Material	White plastic coated galvanised steel	
Shelves	10	
Certification		
Internal Dimension	1100 x 454 x 1325 mm	1100 x 684 x 1325 mm
External Dimension	1220 x 642 x 1885 mm	1220 x 872 x 1885 mm
Weight	158 kg	207 kg
Electrical Requirement	AC220 V/ 50Hz	
Catalogue No.	500300502	500301002



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States