

Pharmacy Refrigerator IUPR-1306



Upright Pharmacy Refrigerator 2 to 8°C IUPR-1306 with transparent door design for easy observation with temperature range 2 to 8°C, and forced air R600a and R134a refrigeration system with air cooler chiller distribution delivers cold air to different levels within each side of the storage chamber, ensuring consistent cabinet wide cooling with adjustable shelves are adjustable to any height and to fulfill different requests.

Ergonomic design offers a clear view of stored items over large framed openings. The thin shape allows for cool recovery of your products, without the concern for swinging door clearances. Customers can choose from two types of interior shelving to suit their needs, security and safety lock to avoid unwanted access, provide storage space for easy sorting of a variety of blood products, and castes design and LED interior light.

Technical Specification:

Model	IUPR-1306
Capacity (liters)	306 L
Cabinet Type (Upright/Chest)	Upright
Temperature Range (°C)	2 to 8 °C
Refrigeration System	R134a
Defrost	Auto
Alarms	High and low temperature alarms, Sensor error alarm, Power failure alarm, Door ajar alarm
Overall Dimensions (WxDxH) cm	640×580×1920 mm
Net Weight (kg)	95 kg
Gross Weight (kg)	110 kg
Electrical Requirements	220V/50Hz
Catalogue Number	500351502

Technical Data:

Model	IUPR-1130	IUPR-1236	IUPR-1306	IUPR-1316	IUPR-1406	IUPR-1416
Capacity	130 L	236 L	306 L	316 L	406 L	416 L
Type	2-8° Pharmacy Refrigerator					
Cabinet Type	Upright					
Temperature Range	2 to 8 °					
Cooling	5°					
Alarms	High and low temperature alarms, Sensor error alarm, Power failure alarm, Door ajar alarm					
Refrigerant	R600a	R134a	R600a	R134a	R600a	R600a
Cooling Type	Forced Air Cooling					
Defrosting	Auto					
Ambient Temp.			10 to 32°		10 to 32°	
Humidity	<60°					
Certification						
Power Consumption	150 W	150 W	220 W	220 W	230 W	190 W
Internal Dimension	418×390×761 mm	500×420×1090 mm	540×425×1380 mm	540×425×1380 mm	685×425×1380 mm	685×429×1380 mm
External Dimension	500×503×1292 mm	580×560×1785 mm	640×580×1920 mm	640×582×1880 mm	780×580×1920 mm	780×580×1920 mm
Net Weight	53 kg	63 kg	95 kg	92.5 kg	115 kg	120 kg
Gross Weight	57 kg	65 kg	110 kg	100 kg	128 kg	133 kg
Electrical Requirement	220V/50Hz					
Catalogue No.	500350502	500351002	500351502	500352002	500352502	500353002



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States