

# **Semi Autoclavable Fixed Volume Micropipette SFG-5000**



Semi Autoclavable Fixed Volume Micropipette SFG-5000 is designed for many different types of routine laboratory work. Build-in ergonomic portable, fits comfortably in user hand.

The basic features of fixed volume pipettes are ideal for clinical diagnostics, control analysis, etc.

The SFG-5000 Fixed Volume Semi Autoclavable Micropipettes is available in capacity of 5000 µl.

## Technical Specification:

|                 |            |
|-----------------|------------|
| Model           | SFG-5000   |
| Volume          | 5000 µl    |
| Accuracy ±%     | 0.6        |
| Accuracy ± (µl) | 30         |
| CV ±%           | 0.2        |
| C V ± (µl)      | 10         |
| Catalogue No.   | 4000570-04 |

## Technical Data:

| Model     | Volume (µl) | Accuracy |      | CV  |      | Catalogue No. |
|-----------|-------------|----------|------|-----|------|---------------|
|           |             | ±%       | ±µl  | ±%  | ±µl  |               |
| SFG-1     | 1           | 5        | 0.05 | 5   | 0.05 | 4000505-04    |
| SFG-2     | 2           | 4        | 0.08 | 2   | 0.04 | 40005010-04   |
| SFG-5     | 5           | 2        | 0.1  | 1   | 0.05 | 4000515-04    |
| SFG-10    | 10          | 1        | 0.1  | 0.5 | 0.05 | 4000520-04    |
| SFG-20    | 20          | 0.8      | 0.16 | 0.4 | 0.08 | 4000525-04    |
| SFG-25    | 25          | 0.8      | 0.2  | 0.4 | 0.1  | 4000530-04    |
| SFG-50    | 50          | 0.8      | 0.4  | 0.4 | 0.2  | 4000535-04    |
| SFG-100   | 100         | 0.6      | 0.6  | 0.2 | 0.2  | 400054004     |
| SFG-200   | 200         | 0.6      | 1.2  | 0.2 | 0.4  | 4000545-04    |
| SFG-250   | 250         | 0.6      | 1.5  | 0.2 | 0.5  | 4000550-04    |
| SFG-500   | 500         | 0.6      | 3    | 0.2 | 1    | 4000555-04    |
| SFG-1000  | 1000        | 0.6      | 6    | 0.2 | 2    | 4000560-04    |
| SFG-2000  | 2000        | 0.6      | 12   | 0.2 | 4    | 4000565-04    |
| SFG-5000  | 5000        | 0.6      | 30   | 0.2 | 10   | 4000570-04    |
| SFG-10000 | 10000       | 0.4      | 40   | 0.2 | 20   | 4000575-04    |



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States