

Semi Autoclavable Fixed Volume Micropipette SFG-10000



Semi Autoclavable Fixed Volume Micropipette SFG-10000 is designed for many different types of routine laboratory work. Build-in ergonomic portable, fits comfortably in user hand.

The basic features of fixed volume pipettes are ideal for clinical diagnostics, control analysis, etc.

The SFG-10000 Fixed Volume Semi Autoclavable Micropipettes is available in capacity of 10000 μ l.

Technical Specification:

Model	SFG-10000
Volume	10000 μ l
Accuracy $\pm\%$	0.6
Accuracy $\pm$ (μ l)	40
CV $\pm\%$	0.2
C V $\pm$ (μ l)	20
Catalogue No.	4000575-04

Technical Data:

Model	Volume (μ l)	Accuracy		CV		Catalogue No.
		$\pm\%$	$\pm\mu$ l	$\pm\%$	$\pm\mu$ l	
SFG-1	1	5	0.05	5	0.05	4000505-04
SFG-2	2	4	0.08	2	0.04	40005010-04
SFG-5	5	2	0.1	1	0.05	4000515-04
SFG-10	10	1	0.1	0.5	0.05	4000520-04
SFG-20	20	0.8	0.16	0.4	0.08	4000525-04
SFG-25	25	0.8	0.2	0.4	0.1	4000530-04
SFG-50	50	0.8	0.4	0.4	0.2	4000535-04
SFG-100	100	0.6	0.6	0.2	0.2	400054004
SFG-200	200	0.6	1.2	0.2	0.4	4000545-04
SFG-250	250	0.6	1.5	0.2	0.5	4000550-04
SFG-500	500	0.6	3	0.2	1	4000555-04
SFG-1000	1000	0.6	6	0.2	2	4000560-04
SFG-2000	2000	0.6	12	0.2	4	4000565-04
SFG-5000	5000	0.6	30	0.2	10	4000570-04
SFG-10000	10000	0.4	40	0.2	20	4000575-04



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States