

Fully Autoclavable Fixed Volume Micropipette AFG-1



Fully Autoclavable Fixed Volume Micropipette AFG-1 is designed for many different types of routine laboratory work. Build-in ergonomic portable, fits comfortably in user hand. The basic features of fixed volume pipettes are ideal for clinical diagnostics, control analysis, etc.

The AFG-1 Fixed Volume "Fully Autoclavable" Micropipettes are available in all capacities from 1 µl.

Technical Specification:

Model	AFG-1
Volume	1 µl
Accuracy ±%	5
Accuracy ± (µl)	0.1
CV ±%	5
C V ± (µl)	0.1
Catalogue No.	4001005-04

Technical Data:

Model	Volume (µl)	Accuracy		CV		Catalogue No.
		±%	±µl	±%	±µl	
AFG-1	1	5	0.1	5	0.1	4001005-04
AFG-2	2	4	0.1	2	0	4001010-04
AFG-5	5	2	0.1	1	0.1	4001015-04
AFG-10	10	1	0.1	1	0.1	4001020-04
AFG-20	20	1	0.2	0	0.1	4001025-04
AFG-25	25	1	0.2	0	0.1	4001030-04
AFG-50	50	1	0.4	0	0.2	4001035-04
AFG-100	100	1	0.6	0	0.2	4001040-04
AFG-200	200	1	1.2	0	0.4	4001045-04
AFG-250	250	1	1.5	0	0.5	4001050-04
AFG-500	500	1	3	0	1	4001055-04
AFG-1000	1000	1	6	0	2	4001060-04
AFG-2000	2000	1	12	0	4	4001065-04
AFG-5000	5000	1	30	0	10	4001070-04
AFG-10000	10000	0	40	0	20	4001075-04



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States