

# Fully Autoclavable Fixed Volume Micropipette AFQ-10



Fully Autoclavable Fixed Volume Micropipette AFQ-10 is fully autoclavable, without disassemble or dismantle. Used to dispense sample measurement between 10 µl volume ranges.

Additional features maintenance-free seals, integral tip ejector for access into narrow tubes: single-channel pipettors, calibrated in accordance ISO/DIS 8655 and each pipettor passes a quality control calibration test and is supplied with a certificate and feature Fully autoclavable, portable ergonomic design, clear volume display, Easy and simple calibration, avoid repetitive strain, provide individual test certificate ISO 8655/DIN 12650 with a limited 3-year warranty.

## Technical Specification:

|               |            |
|---------------|------------|
| Model         | AFQ-10     |
| Volume        | 10 µl      |
| Precision     | 0.8%       |
| Accuracy ±%   | 0.8%       |
| Catalogue No. | 4001510-04 |

## Technical Data:

| Model    | Volume (µl) | Test Volume | Inaccuracy | Imprecision | Catalogue No. |
|----------|-------------|-------------|------------|-------------|---------------|
| AFQ-5    | 5µl         | 5µl         | 1.3%       | 1.2%        | 4001505-04    |
| AFQ-10   | 10µl        | 10µl        | 0.8%       | 0.8%        | 4001510-04    |
| AFQ-20   | 20µl        | 20µl        | 0.6%       | 0.5%        | 4001515-04    |
| AFQ-25   | 250µl       | 250µl       | 0.5%       | 0.3%        | 4001520-04    |
| AFQ-50   | 50µl        | 50µl        | 0.5%       | 0.3%        | 4001525-04    |
| AFQ-100  | 100µl       | 100µl       | 0.5%       | 0.3%        | 4001530-04    |
| AFQ-200  | 200µl       | 200µl       | 0.4% 0     | .2%         | 4001535-04    |
| AFQ-250  | 250µl       | 250µl       | 0.4%       | 0.2%        | 4001540-04    |
| AFQ-500  | 500µl       | 500µl       | 0.3%       | 0.2%        | 4001545-04    |
| AFQ-1000 | 1000µl      | 1000µl      | 0.3%       | 0.2%        | 4001550-04    |
| AFQ-2000 | 2000µl      | 2000µl      | 0.3%       | 0.2%        | 4001555-04    |
| AFQ-5000 | 5000µl      | 5000µl      | 0.3%       | 0.15%       | 4001560-04    |



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States