

# Universal Trays 755 x 481 mm



The Universal Trays employ a drilled grid for individual loadings of test tube holders, microtitre plates, and deep well plates.

Suitable for	OS-3050	OS-3100	OS-3150	OS-3200
Dimensions (W x D in mm)	320 x 260	350 x 350	450 x 450	755 x 481
<b>Max. Load of Beaker Holders (kg)</b>				
Beaker Volume				
50 ml	21	36	55	99
100 ml	14	28	44	74
250 ml	8	13	24	39
300 ml	8	13	24	39
500 ml	5	10	16	25
1 L	3	5	8	14
2 L	1	4	5	9
2.8 L	1	2	4	6
4 L	-	1	2	4
<b>Max. Load of Separatory Funnel Holders (kg)</b>				
Separatory Funnel Volume				
250 ml	4	4	6	11
500 ml	3	3	4	9
1 L	2	2	3	5
2 L	-	1	2	4



**Iris Analytical Ltd.**

66 W Flagler St #900 Miami, Florida 33130 United States