

Multi-Purpose Tray for Orbital Shakers OS-3050



The multi-purpose tray is made of stainless steel and can be used for a diverse range of unusual specifications.

Suitable for	OS-3050	OS-3100	OS-3150	OS-3200
Dimensions (W x D in mm)	320 x 260	350 x 350	450 x 450	755 x 481
Max. Load of Beaker Holders (kg)				
Beaker Volume				
50 ml	12	16	25	40
100 ml	6	9	16	28
250 ml	3	4	9	15
300 ml	3	4	9	15
500 ml	3	4	5	10
1 L	2	2	4	8
2 L	-	1	2	4
2.8 L	-	1	2	4



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States