

# Multi-Purpose Tray for Orbital Shakers OS-3150



The multi-purpose tray is made of stainless steel and can be used for a diverse range of unusual specifications.

| Suitable for                     | OS-3050   | OS-3100   | OS-3150   | OS-3200   |
|----------------------------------|-----------|-----------|-----------|-----------|
| Dimensions (W x D in mm)         | 320 x 260 | 350 x 350 | 450 x 450 | 755 x 481 |
| Max. Load of Beaker Holders (kg) |           |           |           |           |
| Beaker Volume                    |           |           |           |           |
| 50 ml                            | 12        | 16        | 25        | 40        |
| 100 ml                           | 6         | 9         | 16        | 28        |
| 250 ml                           | 3         | 4         | 9         | 15        |
| 300 ml                           | 3         | 4         | 9         | 15        |
| 500 ml                           | 3         | 4         | 5         | 10        |
| 1 L                              | 2         | 2         | 4         | 8         |
| 2 L                              | -         | 1         | 2         | 4         |
| 2.8 L                            | -         | 1         | 2         | 4         |



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States