

Unstirred Water Bath WB-10



- Tank is using seamless stainless steel and outside body is powder painting
- Screen-touch panel, micro-computer PID digital display
- Overheat protector
- The water bath with circulator using electric heat tube to unify temperature inside
- The sensor is PT-100Ω
- Gable cover (AC-100) sold separately for your option.
- Heater of stainless steel

Technical Data:

Model	WS-100		WS-200	
Temperature Range (?)	Ambient+5?~100?			
Temp. Accuracy (?)	±0.1?			
Temp. Control	PIDControl			
Temp. Display	LEDDisplay			
Volume (liters)	10 L		20 L	
Inside Dimensions (WxDxH) mm	295 x 235 x 150 mm		495 x 295 x 150 mm	
Outside Dimensions (WxDxH) mm	340 x 280 x 315 mm		540 x 340 x 315 mm	
Net Weight (kg)	8kg		12kg	
Gross Weight (kg)				
Electrical Requirements	110V/60Hz	220V/50Hz	110V/60Hz	220V/50Hz
Catalogue No.				



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States