

Water Bath WS-100



- Tank is using seamless stainless steel and outside body is powder painting
- Screen-touch panel, micro-computer PID digital display
- Overheat protector
- The water bath with circulator using electric heat tube to unify temperature inside
- Add exterior circulator as require
- The sensor is PT-100Ω
- Gable cover (AC-100) sold separately for your option.
- Heater of stainless steel
- External Circulating Pump is Available
- Heating circulator sold separately for your option
- Optional Water level protector
- Optional Timer

Technical Data:

Model	WS-100		WS-200	
Temp. Range (?)	Ambient +5?~100?			
Temp. Accuracy (?)	±0.1?			
Temp. Control	PID Control			
Temp. Display	LED			
Volume	10 L		20 L	
Inside Dimension (WxDxH) mm	295 x 235 x 150		495 x 295 x 150	
Outside Dimension (WxDxH) mm	340 x 280 x 310		540 x 340 x 310	
Net Weight (kg)	7 kg		9 kg	
Gross Weight (kg)				
Electrical Requirements	110V / 60Hz	220V / 50Hz	110V / 60Hz	220V / 50Hz
Catalogue No.				



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States