

Wide Field Zoom Stereo Microscope IMS-3260



Wide Field Zoom Stereo Microscope IMS-3260 is designed for large viewing field without eyepiece for sharp and stereoscopic images from the lens in microscope head and has excellent objective to choose 2X, 4X, 6X, with optional choice 8X, 10X, 15X and 6XSL objectives. IMS-3260 model is capable of providing strong images with a broad field of view, the multipurpose and perfect for basic studies of samples properties. Come with LED lights, allowing samples to be illuminated for a more comprehensive field of view. Suitable for Industrial, Medical and Scientific Research Areas, circuit board repair, circuit board inspection, surface mount technology work, electronics inspection, coin collecting, gemology and gemstone setting, engraving, and repair and inspection of small parts. Provide upright, un-reversed 3 D images. Used to examine the huge specimen under 4X-12X zoom magnification objective.

Technical Specification:

Model	IMS-3260
Type	Wide Field Zoom Stereo Microscope
Viewing Head	2× high quality optical head without eyepieces
Focusing	Coarse focus with focus knob on stand
Objective	2X208mm, FOV: 68mm, 4X98mm, FOV: 34mm, 6X80mm, FOV: 22.7mm
Light Source	3W built-in LED illumination with external control box
Magnification	4X-12X
Pillar Stand	Column height 580mm, moving range 330mm, Base 400×300×20mm, stage size: 210×190mm
Catalogue No.	440205-10-04

Technical Data:

Model	IMS-3255	IMS-3260	IMS-3265
Type	Wide Field Zoom Stereo Microscope		
Viewing Head	2× high quality optical head without eyepieces		
Magnification	4X-12X		
Focusing	Coarse focus with focus knob on stand		Focus with handle on head
Eyepiece			
Objective	2X208mm, FOV: 68mm, 4X98mm, FOV: 34mm, 6X80mm, FOV: 22.7mm, (optional: 8X 58mm, FOV: 17mm, 10X46mm, FOV: 13.6mm, 15X50mm, FOV: 9.1mm, 6XSL(Super Long working distance), 115mm, FOV: 22.7mm)		
Pillar Stand	Column height 500mm, moving range 175mm, Base 400×300×20mm	Column height 580mm, moving range 330mm, Base 400×300×20mm, stage size: 210×190mm	Universal stand can be fixed on table
Light Source	3W built-in LED illumination with external control box		
Certification			
Packing Dimension			
Packing Weight			
Electrical Requirement			
Catalogue No.	440205-05-04	440205-10-04	440205-15-04



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States