

Analytical Balances 0.1mg IBA-2100(110V)



Analytical Balances 0.1mg IBA-2100 feature weighing capacities of 100 g, and readabilities of 0.1 mg. Our Analytical Balances are designed with for pharmaceutical, laboratory and regulated weighing in a variety of cell and tissue culture facilities that require the highest accuracy in analytical weighing processes.

Iris Analytical Balances are designed for routine weighing applications in your workplace, offering accuracy and repeatability during essential weighing applications in education, laboratory and, industrial sector. Iris analytical balances deliver competitive performance at an economical price, featuring RS232 connectivity for easy communication, and a backlit display and a simple interface for uncomplicated operation, these balances are perfectly designed for your workplace.

Standard features include Basic functionality for routine weighing applications. Iris analytical balances are equipped with three essential weighing modes, RS232 connectivity for data transfer and storage and internal calibration (in selected models) making them ideal for routine weighing applications.

Designed for uncomplicated operation with easy-to-use display and interface Equipped with an easy-to-read, bright backlit display and a simple user interface, they are incredibly easy to operate, with almost no training required.

Smart design and durable construction Iris analytical balances have a small footprint which saves on benchtop space while providing a large weighing surface. They are durably constructed, and feature a stainless-steel pan to withstand day-to-day use in the workplace.

Technical Specification:

Model	IBA-2100-110V
Capacity	0 to 100 g
Repeatability	0.1 mg
Resolution	±0.2 mg
Display	LCD
Pan Dimensions	80 mm
Construction Material	Base: cast Al Alloy STAINLESS: platter
Linearity	±0.2 mg
Calibration	external or internal
Measuring Parameters	g/mg/ct/oz
Interface	RS232, USB
Dimensions	310x 212x350 mm
Net Weight	50 kg
Electrical Requirements	AC110V
Catalogue No.	680250501

Technical Data:

Model	IBA-2100	IBA-2200	IBA-2300	IBA-2400	IBA-2500	IBA-2600	IBA-2700	IBA-2800	IBA-2900	IBA-2980	
Capacity	0 to 100 g	0 to 120 g	0 to 160 g	0 to 180 g	0 to 200 g	0 to 210 g	0 to 220 g	0 to 300 g	0 to 310 g	0 to 320 g	
Resolution	±0.2 mg										
Display	LCD										
Resolution	0.1 mg										
Pan Dimension	80 mm										
Material	Base: cast Al Alloy STAINLESS: platter										
Linearity	±0.2 mg										
Calibration	external or internal										
Measuring Parameter	g/mg/ct/oz										
Interface	RS232, USB										
Dimension	310x 212x350 mm										
Weight	50 kg										
Electrical Requirement	AC110V	240V	NA								
Catalogue No.	680250501	680250502	6802510-04	6802515-04	6802520-04	6802525-04	6802530-04	6802535-04	6802540-04	6802545-04	6802550-04



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States