LCD Digital Rotary Evaporator is an essential instrument in organic and pharmaceutical labs to remove low boiling solvents from a flask and to get the requested compound out of the solution by heating a rotating flask that contains the user sample in a reduced pressure environment.

- Large digital LCD screen displays heating temperature, rotation speed and time
- Motor lift configurations
- Wide temperature range of the 5L heating bath
- Speed range from 20 to 280 rpm, and interval operation in clockwise and counter-clockwise directions
- Patented condenser (cooling surface 1,500 cm2) with excellent cooling effect
- Easy evaporating flask removal thanks to the ejector
- Chemical-resistant double PTFE system and patented pressure spring provide excellent sealing
- PC connectivity via USB interface

Rotary Vacuum Evaporators

Rotary Vacuum Evaporators





IRIS Analytical Inc.

cs@irisanalytical.com +1 706 760 4500 www.irisanalytical.com