IRIS ANALYTICAL



MICROSCOPES

cs@irisanalytical.com


+1 706 760 4500

www.irisanalytical.com

Stereo Microscope IMS-1300 Series suitable for school education, circuit board inspection, and repair, surface mount technology work, electronics inspection, coin collecting, gemology and gemstone setting, engraving, inspection and repair of small parts. The binocular viewing head and V-shaped objective. The stand has the same ergonomic design with the incident lighting systems are rheostat controlled by 3W LED lamp that allows the illumination to be adjusted to suit the sample type.

- Long Working Distance: 1X,230 mm objective to create enough space for the uses
- Extra wild field WF10X20 mm eyepieces
- ELWD Condenser NA 0.3, LWD 72 mm
- Incident 3W LED light source with incident illumination
- Different filters including Blue, Green and Ground glass filters
- Binocular 45° inclined and a 360° rotatable viewing head with diopter adjustment ±5dp on left ocular-tube and tube diameter 23.2 mm
- Small foot print and light weight
- Easy operating with ergonomic design
- Optional accessories Eyepiece, objective

IRIS ANALYTICAL



IRIS Analytical Inc.

cs@irisanalytical.com +1 706 760 4500 www.irisanalytical.com