IRIS ANALYTICAL





STEAM STERILIZERS/AUTOCLAVES

cs@irisanalytical.com +1 706 760 4500 www.irisanalytical.com


Tabletop Autoclave ITA-400 Series

The ITA-400 front loading ‘Benchtop’ autoclave is pulse vacuum type and range is available in single sizes of 24 litres chamber capacity. These units require no plumbing, and only a standard single phase power supply. Build in 7 Sterilization cycles + 2 User definable cycles + 2 Test programs.

- LCD display, 7 Sterilization cycles + 2 User definable cycles + 2 Test programs
- Full-Automatic operation, check and monitor by computer
- printer to note the sterilization process in details
- Automatic safety door lock system
- Pulse Vacuum
- The water setting is fully automatic
- Connector for Printer sets, SD Cardas an optional accessories
- PID control

Model	ITA-424
Volume of Chamber	24 L
Air Expulsion	Pulse Vacuum
Sterilization Temp	121? (1.1 kg/cm²) / 134? (2.1 kg/cm²)
Dry Function	Adjustable 3 to 30 minutes
Water Reservoir	Supply tank : 4.5L / Discharge tank : 2.5L
Sterilization Programs	7 Sterilization cycles + 2 User definable cycles + 2 Test programs
Standard Accessory	Sterilization frame, sterilization network disk (375x220x20mm x 4 drive, Note: 3Kg / plate) SUS304 sterilization clip x1 branch network disk removed
Optional Accessory	Printer sets, SD Card
Power Consumption	2000 W
Net Weight	57 kg
Gross Weight	60 kg
Size of Chamber	260x455 mm
External Dimensions	560x665x400 mm
Electrical Requirement	AC 230V±10% 50HZ / 60HZ, 9A
Catalogue No.	580150502

IRIS ANALYTICAL



IRIS Analytical Inc.

cs@irisanalytical.com +1 706 760 4500 www.irisanalytical.com