IRIS ANALYTICAL

CENTRIFUGES



cs@irisanalytical.com

+1 706 760 4500


www.irisanalytical.com

Portable Centrifuge IPC-400 Series is perfect for quick spin downs the Multi-Rotor can be placed 12 microtubes and /or 4(8X0.2ml) PCR strips and the Strip-Rotor can be placed 4 (12X0.2ml) PCR strips. The Max speed can reach 12000rpm/min. The IPC-400 Series is a personal centrifuge with unmatched capacity and flexibility.

- Time from 1sec-999min
- Multi-Rotor for tubes and Strips
- LCD display of time and speed
- The double capacity of minicentrifuge
- Adjustable speed from 1000rpm-12000rpm

Model	IPC-410		IPC-420	
Max. capacity	12X1.5ml / 2.0ml tubes & 4xPCR strips(8X0.2ml) ; 4xPCR strips(12X0.2ml); 4X5ml;		12X1.5ml / 2.0ml tubes & 4xPCR strips(8X0.2ml) ; 4xPCR strips(12X0.2ml); 4X5ml;	
Max. speed rpm (Angle rotor)	7,000 rpm		12,000 rpm	
Speed increment	1000 rpm		1000 rpm	
Max. RCF[×g]	3240 x g		9151 x g	
Rotor	24 Capillary haematocrit Rotor		24 Capillary haematocrit Rotor	
Display	LCD display			
Time Range	1sec-999min			
Electrical Requirements	110V/60Hz	220V/50Hz	110V/506z	110V/60Hz
Catalogue No.	140350501	140350502	140351001	140351001

IRIS ANALYTICAL



IRIS Analytical Inc.

cs@irisanalytical.com +1 706 760 4500 www.irisanalytical.com