IRIS ANALYTICAL

CENTRIFUGES



cs@irisanalytical.com

+1 706 760 4500

www.irisanalytical.com


Blood type ID Card Centrifuge IIDC-2000blood card centrifuge especially used for blood group testing, microcolumn gel immune detection. Spins sample volumes up to12 and 24 cards micro gel rotor with maximum speed 3840 rpm. Driven by a brushless motor. Auto-balancing centrifuge for blood bank is particularly designed for the quick centrifuge in blood bank with storages up to 9 program records for users. For safety, we introduced an electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function.

- Specifically for blood group test, micro-column gel immune detection
- The machine adopts brushless motor driving,micro-computer control
- Electric lid lock, super speed, and self-test protection
- Low noise, without carbon dust pollution
- Professional program design Can be run directly Without tedious parameter settings
- Maximum capacity of up to 12 and 24 cards micro gel rotor

Blood group card dedicated program automatically

Timer Range	Speed	RCF
0 to 2 min	780 rpm	75xg
2 to 3 min	1280 rpm	200xg
3 to 10 min	3840 rpm	1790xg

IRIS ANALYTICAL



IRIS Analytical Inc.

cs@irisanalytical.com +1 706 760 4500 www.irisanalytical.com