IRIS ANALYTICAL

CENTRIFUGES



cs@irisanalytical.com


+1 706 760 4500

www.irisanalytical.com

Blood type ID Card Centrifuge IIDC-3000 with standard programs are built-in for blood type determination and all kinds of experiments of serology, which standardize the test and research of clinical. Users can set the programs with a simple touch screen button. Suitable for blood group testing, microcolumn gel immune detection. Spins sample volumes up to 24 cards micro gel rotor with a maximum speed of 3000 rpm. Driven by a DC brushless motor. Auto-balancing centrifuge for blood bank is particularly designed for the quick centrifuge in blood bank with storages up to 9 program records for users. For safety, we introduced an electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function.

- Specifically for blood group test, micro-column gel immune detection
- The machine adopts DC brushless motor driving,micro-computer control
- Electric lid lock, super speed, and self-test protection
- Low noise, without carbon dust pollution
- Professional program design Can be run directly Without tedious parameter settings
- Maximum capacity of up to 24 cards micro gel rotor

IRIS ANALYTICAL



IRIS Analytical Inc.

cs@irisanalytical.com +1 706 760 4500 www.irisanalytical.com