IRIS ANALYTICAL

REFRIGERATORS & ULTRA-LOW FREEZERS



cs@irisanalytical.com


+1 706 760 4500

www.irisanalytical.com



Ultra-Low Upright Freezers -86℃ IUUF-2000 Series design with a cascade refrigeration system that uses two kinds of refrigerants having different boiling points, which run through their own independent freezing cycle and are joined by a heat exchanger. Features with filter for easy cleaning and single hand openable, Gasket heating, leveling caster, etc. Available with maximum capacity up to 390 L490 L, 650 L with a temperature range of -50 to-86°C.

IRIS ANALYTICAL



IRIS Analytical Inc.

cs@irisanalytical.com +1 706 760 4500 www.irisanalytical.com