IRIS ANALYTICAL

LAMINAR FLOW CABINETS




cs@irisanalytical.com

+1 706 760 4500

www.irisanalytical.com

Vertical Laminar Flow Cabinets LVC-1000 Series provide a particle-free working environment by taking air through a filtration system and exhausting it across a work surface in a laminar or unidirectional air stream. The laminar flow cabinet is enclosed on the sides and kept under constant positive pressure in order to prevent the infiltration of contaminated room air. Laminar flow cabinet is widely used in medical research laboratories, hospitals, manufacturing facilities, and other research and production environments.

IRIS ANALYTICAL



IRIS Analytical Inc.

cs@irisanalytical.com +1 706 760 4500 www.irisanalytical.com