

UV-VIS SPECTROPHOTOMETERS



cs@irisanalytical.com


+1 706 760 4500

www.irisanalytical.com

Single Beam UV-VIS Spectrophotometer IUS-1000 offers superb application functions with accuracy and precision. It is the best instrument for various applications. Provide 1200 line grating single beam to ensure high resolution with wavelength availability of 190 to 1000 nm and 4 nm bandwidth, tungsten, and D2 lamp. A useful function is an absorption, transmittance, and wavelength with automatic switching mode.

- Equipped with a microprocessor to set parameter
- 1200 line single beam grating which ensures high resolution
- Designed with 128×64-inch and 320×240-inch LCD display
- Variable spectral bandwidth 0.5, 1, 1.8, 2, 4 and 5 nm
- Die-cast aluminum base materials
- Single beam spectrophotometer with highly stable optics
- Low stray light with ≤ 0.2% T
- Automatic setting mode
- Designed with multi-lamp source
- Offer high-class grating hermetic light path to ensure low stray light
- Support USB Port & Parallel Port (Printer)

IRIS ANALYTICAL



IRIS Analytical Inc.

cs@irisanalytical.com +1 706 760 4500 www.irisanalytical.com