IRIS ANALYTICAL

Filtration **Products**



cs@irisanalytical.com


+1 706 760 4500

www.irisanalytical.com

Vacuum Filtration Manifold IFM-300 Series provides independent control in every branch for the vacuum and venting process. Connect 3 and 6 filters vessels to one vacuum source. This contains flasks, filter cup adaptors, and vacuum and autoclavable manifold and provides filters for many applications. The filtration manifold usages a patented spin-lock assembly that enables quick and stable installation without a clamp. It has been used for all kinds of liquid filtration, biological laboratories, chemical industry, microbiological detection.

- The stainless steel funnels that can be flame sterilized
- Build with patented spin-lock design
- Quick, easy and safe to service, fast-assembly and disassembly
- Requires electricity for operating
- Large capacity in small footprint, exchangeable funnel adaptor design permits user to fit their desired funnel simply
- Cost-effective high capacity system
- Available capacity 100, 300 and 500 ml
- Accessories: stainless steel manifold, 3 x 100/300/500 ml Stainless Steel Filter Cups and 6 x 100/300/500 ml Stainless Steel Filter Cups

IRIS ANALYTICAL



IRIS Analytical Inc.

cs@irisanalytical.com +1 706 760 4500 www.irisanalytical.com