IRIS ANALYTICAL





WATER/ WASTE WATER TESTING

cs@irisanalytical.com


+1 706 760 4500

www.irisanalytical.com

Portable Ion Meter IM-200is a multi-used, strong, and accurate instrument. The easy features operation menu with easy to understand prompts and electrode/standard error checking. These meters feature many advanced options, such as programmable stability measures, programmable standardization reading delay times, multi-channel operation, fast reading update rates of twice per second for all channels, programmable alarms, storage programmable data up to 500 to 1400 data points and RS232, USB serial interface for controlling the meter and obtaining data. The meter operating code can easily be upgraded as new features are made available.

- Designed with 3 reading modes and read function
- Easily read the sample changing in the measurement process
- Power-off safety
- Designed with Rubber foot
- Simple calibration up to 5 point calibration with automatic
- Functional measuring parameters including pX/Ion concentration/mV/ Temperature
- Large, well-structured TFT Touch screen and LCD display with navigation
- CE Certified
- Manual temperature compensation and temperature coefficient
- Maximum storage: 1000 to 1400 sets
- Auto switch off mode
- Advanced functions in the identification of standard buffer and ID management
- Compiled with micro USB and GLP
- Li-on battery and AA battery
- Available standard accessories: Fluoride ISE, Electrode,TATC Probe, Stirrer

IRIS ANALYTICAL



IRIS Analytical Inc.

cs@irisanalytical.com +1 706 760 4500 www.irisanalytical.com