

Homogenizers & Dispersers


IRIS ANALYTICAL

cs@irisanalytical.com
+1 706 760 4500
www.irisanalytical.com

Overhead Stirrers OS-1000 Series suitable for all kinds of conventional mixing applications. The simple and intuitive operation can significantly improve the work efficiency in the laboratory. rent optional impellers.

- Directly connectable to the computer via RS-232 cable for remote control
- Equipped with digital input/output ports and foot pedal port for various applications
- Suitable for stirring in reaction systems or stirring applications linked to other equipment
- Rotation speed can be lowered to 20 rpm for very sensitive stirring processes
- Brushless DC motor, which is very suitable for long-term experimental applications
- Stable and accurate stirring process due to latest micro-processor technology
- Totally enclosed and compact casing
- Stable and quiet working process
- Internal overload protection
- Adjustable impeller shaft for different heights
- Soft start at low speed guarantees a smooth and safe stirring process
- Suitable for applications in various environments
- Comes with a reversible rotation function
- The standard version package entails the overhead stirrer, a stand, rod, and clamp, as well as a stainless steel impeller
- Availability of a wide range of different

IRIS ANALYTICAL



IRIS Analytical Inc.

cs@irisanalytical.com +1 706 760 4500 www.irisanalytical.com