IRIS ANALYTICAL



DRY BLOCK HEATER/THERMO MIXERS

cs@irisanalytical.com | +1 706 760 4500 | www.irisanalytical.com


The Digital Dry Bath Incubator DI-300 series are ideal for incubation and activation of cultures, enzyme reactions, immunoassays, melting/boiling points and a wide variety of other laboratory procedures. DI-300DE is an innovative design accommodating two blocks with independent digital temperature control. Each block can be set at different temperatures - ideal for multiple users or for applications where samples have to be transferred between two temperatures very quickly.

Features:

- Exceptional temperature uniformity and accuracy
- Wide range of interchangeable aluminium blocks
- Heat up to 105℃, optional external temperature probe
- Blocks as accessories for all applications -- tubes, vials and micro plates
- Block extraction tool is supplied, allowing blocks to be removed easily
- Auto-stop and sound-alert after program completed

Model	DI-300		DI-300D		DI-300DE	
Temperature Range	RT+5°C to 105°C		-RT+5°C to 105°C		RT+5°C to 105°C	
Timing Range	1min to 99h59min		1min to 99h59min		1min to 99h59min	
Temperature Stability@100°C	?0.5°C		?0.5°C		?0.5°C	
Uniformity: within the block@37°C	±0.2°C		±0.2°C		±0.2°C	
Across similar block@37°C	----		±0.3°C		----	
Heating Time(25 to100°C)	?15min		?15min		?15min	
Display Accuracy	±0.1°C		±0.1°C		±0.1°C	
Block Quantity	1		2		2	
Net Weight	2.5kg		4.5kg		4.5kg	
Dimension	250 x 190 x130 mm		365 x 210 x150 mm		365 x 210 x150 mm	
Electrical Requirement	110V/60Hz	220V/50Hz	110V/60Hz	220V/50Hz	110V/60Hz	220V/50Hz
Catalogue No.	1601505-01	1601505-02	1601510-01	1601510-02	1601515-01	1601515-01

IRIS ANALYTICAL



IRIS Analytical Inc.

cs@irisanalytical.com +1 706 760 4500 www.irisanalytical.com