IRIS ANALYTICAL

LAMINAR FLOW CABINETS



cs@irisanalytical.com


+1 706 760 4500

www.irisanalytical.com


Vertical Laminar Flow Cabinet LVC-2000 Series provide a particle-free working environment by taking air through a filtration system and exhausting it across a work surface in a laminar or unidirectional air stream. The laminar flow cabinet is enclosed on the sides and kept under constant positive pressure in order to prevent the infiltration of contaminated room air.

Laminar flow cabinet is widely used in medical research laboratories, hospitals, manufacturing facilities and other research and production environments.

LED Display
Microprocessor control system, LED display.




UV Lamp
Emission of 253.7 nanometers for most efficient decontamination.



Double Sides
Double sides work bench, operators can sit face to face.



IRIS ANALYTICAL



IRIS Analytical Inc.

cs@irisanalytical.com +1 706 760 4500 www.irisanalytical.com