

# PUMPS & ASPIRATORS

cs@irisanalytical.com


+1 706 760 4500

www.irisanalytical.com

**Basic Peristaltic Pump ISPP-1250S** has 0.000829 to 2280 ml/min flow rates and excellent dosing that used for timer function of simple dispensing function along with RS232 standard interface along with optional RS485 Modbus interface communicator.

- Support RS232 and RS485, the pump can be connected with PC to reach dispensing and filling function
- 3 Digital LED display the current speed
- ABS plastic constructing material
- A quick exchange of pump heads without the use of tools
- The strong, powerful, regulated drive unit
- Allow stacking of 4 pumps head for same drive, creating multichannel pump system
- Easy installation and low maintenance condition
- Materials have good rigidity and high precision
- Mechanical keypad to control the start/stop, speed, and direction

IRIS ANALYTICAL




IRIS Analytical Inc.

cs@irisanalytical.com +1 706 760 4500 www.irisanalytical.com