

PUMPS & ASPIRATORS

cs@irisanalytical.com


+1 706 760 4500

www.irisanalytical.com

Compact Peristaltic Pump ICPP-1445-1PH has 0.004 to 149.23 ml/min flow rates and excellent dosing that used for timer function of simple dispensing function along with external control0 to 5V, 4 to 20 mA, and 0 to 10V (OPTIONAL) USB connector, RS485 interface.

- Low power consumption
- LED display the current speed
- ABS plastics constructing material
- A quick exchange of pump heads without the use of tools

IRIS ANALYTICAL



IRIS Analytical Inc.

cs@irisanalytical.com +1 706 760 4500 www.irisanalytical.com