Laboratory Glassware Washers IGW-1200(380V/50Hz)



Laboratory Glassware Washers IGW-1200 suitable for washing and drying narrow-neck with HEPA filter, general purpose glassware washer and has specialized features for contamination-sensitive research especially for laboratory and medical applications.

The IGW-1200 is decontaminating equipment that worked on moist heat to remove contamination and bacteria with a large wash cabinet and short programme cycles. The key features include stainless steel interior, exterior, racks and inserts withstand heavy use and provide protection against corrosion resistance materials and offer internal components peristaltic pumps and seals laboratory-grade for excellent durability of instruments.

Technical Specification:

Model	IGW-1200
Functions	Washing
Capacity	176 L
Circulating System	0 to 400 L/min
Materials	Internal Panel: 316L, External Panel: 304L
Peristaltic pump	2 to 4 pcs
Programs	8 Standard, 2 user defined
Interface	RS232, SD card, printer, pure booster pump
Heating Power	4 KW/9 KW
Dimensions	950×605×750 mm
Electrical Requirements	380V/50Hz
Catalogue No.	370050503

Technical Data:

Model	IGW-1200		IGW-1400	
Difference	Washing		Washing + Drying with HEPA filter	
Capacity	176 L		176 L	
Circulating System	0 to 400 L/min		0 to 400 L/min	
Materials	Internal Panel: 316L, External Panel: 304L		Internal Panel: 316L, External Panel: 304L	
Peristaltic pump	2 to 4 pcs		2 to 4 pcs	
Programs	8 Standard, 2 user defined		8 Standard, 2 user defined	
Function	Fault diagnose, electronic lock, automatic power off		Fault diagnose, electronic lock, automatic power off	
Interface	RS232, SD card, printer, pure booster pump		RS232, SD card, printer, pure booster pump	
Heating Power	4 KW/9 KW		4 KW/9 KW	
Dimension	950×605×750 mm		950×925×750 mm	
Electrical Requirement	220V/50Hz	380V/50Hz	220V/50Hz	380V/50Hz
Catalogue No.	370050502	370050503	370051002	370051002



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States