# Semi Autoclavable Fixed Volume Micropipette SFG-25



Semi Autoclavable Fixed Volume Micropipette SFG-25 is designed for many different types of routine laboratory work. Build-in ergonomic portable, fits comfortably in user hand.

The basic features of fixed volume pipettes are ideal for clinical diagnostics, control analysis, etc.

The SFG-25 Fixed Volume Semi Autoclavable Micropipettes is available in capacity of 25 µl.

## Technical Specification:

| | |
|---|---|
| Model | SFG-25 |
| Volume | 25 µl |
| Accuracy ±% | 0.8 |
| Accuracy ± (µl) | 0.2 |
| CV ±% | 0.4 |
| C V ± (µl) | 0.01 |
| Catalogue No. | 4000530-04 |

## Technical Data:

| Model | Volume (µl) | Accuracy | | CV | | Catalogue No. |
|---|---|---|---|---|---|---|
| | | ±% | ±µl | ±% | ±µl | |
| SFG-1 | 1 | 5 | 0.05 | 5 | 0.05 | 4000505-04 |
| SFG-2 | 2 | 4 | 0.08 | 2 | 0.04 | 40005010-04 |
| SFG-5 | 5 | 2 | 0.1 | 1 | 0.05 | 4000515-04 |
| SFG-10 | 10 | 1 | 0.1 | 0.5 | 0.05 | 4000520-04 |
| SFG-20 | 20 | 0.8 | 0.16 | 0.4 | 0.08 | 4000525-04 |
| SFG-25 | 25 | 0.8 | 0.2 | 0.4 | 0.1 | 4000530-04 |
| SFG-50 | 50 | 0.8 | 0.4 | 0.4 | 0.2 | 4000535-04 |
| SFG-100 | 100 | 0.6 | 0.6 | 0.2 | 0.2 | 400054004 |
| SFG-200 | 200 | 0.6 | 1.2 | 0.2 | 0.4 | 4000545-04 |
| SFG-250 | 250 | 0.6 | 1.5 | 0.2 | 0.5 | 4000550-04 |
| SFG-500 | 500 | 0.6 | 3 | 0.2 | 1 | 4000555-04 |
| SFG-1000 | 1000 | 0.6 | 6 | 0.2 | 2 | 4000560-04 |
| SFG-2000 | 2000 | 0.6 | 12 | 0.2 | 4 | 4000565-04 |
| SFG-5000 | 5000 | 0.6 | 30 | 0.2 | 10 | 4000570-04 |
| SFG-10000 | 10000 | 0.4 | 40 | 0.2 | 20 | 4000575-04 |




# Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States