
Multi-Pore Holders - Multi-Layer Holder for OS-3000 Series Orbital shakers



The multi-pore holders can be installed on the sticky plates or the universal trays. All holders can be installed on small pores and deeper holes In total, four different multi-pore holders are available:

1. Multi-Layer Holder
2. Single-Layer Single-Position Holder
3. Single-Layer Three-Position Holder
4. Single-Layer Four-Position Holder

Max. No. of Mountable Holders	Multi-Layer	Four-Position	Three-Position	Single-Position
OS-3050	4	-	1	6
OS-3100	6	-	2	8
OS-3150	10	2	3	15
OS-3200	16	3	5	25



Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States