Conical Centrifuge Tubes, Polypropylene, 50ml



Conical screw cap centrifuge tubes have a large white labeling area and easy-to-read graduations. HDPE* cap has a sealing ring for a leak-proof seal and is flat for easy labeling. Copolymer offers excellent clarity and strength and has the same chemical resistance as polypropylene. Withstands 10,000 x *g* and has a temperature range of -40°C to +120°C. 500 tubes per case.

Iris Analytical Ltd.

66 W Flagler St #900 Miami, Florida 33130 United States